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The screenshot shows the CMDfact software interface. The main window is titled "CMDfact Patientenverwaltung". It is divided into two main sections: "Patientinnen und Patienten" on the left and "Untersuchungen" on the right.

**Patientinnen und Patienten:** This section features a search bar labeled "Suche" with the text "ma" entered. To the right of the search bar is a "Neu einlesen" button with a refresh icon. Below the search bar is a table with two columns: "Name" and "Pat.-Nr.". The table contains two entries: "Mustermann, Martina" with ID "123456" and "Musterfrau, Mara" with ID "133578". A red circle with the number "1" highlights the name "Musterfrau, Mara". A red circle with the number "2" highlights the search bar. Below the table is a button labeled "Diesen Patienten auswählen".

**Untersuchungen:** This section shows a search for ID "123456 |". Below this is a table with three columns: "Datum", "Untersuchung", and "Nr.". The table contains five entries, all dated "21.06.2017": "CMDcheck" (Nr. 1), "CMDstatus" (Nr. 2), "CMDmanu" (Nr. 3), "CMDpain" (Nr. 4), and "CMD DiagnosePilot" (Nr. 5). Below the table are two buttons: "Neuer Patient" and "Patient Bearbeiten". A red circle with the number "3" highlights the "Patient Bearbeiten" button.

At the bottom of the window, there is a status bar with the text "Patientin/Patient Pat.-Nr. - - -" on the left and "Praxis CMD-Centrum Hamburg-Eppendorf" on the right.

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The screenshot shows the 'CMDfact Patientenverwaltung' window. It is divided into two main sections: 'Patientinnen und Patienten' on the left and 'Untersuchungen' on the right. In the 'Patientinnen und Patienten' section, there is a search field labeled 'Suche' with the placeholder text 'Pat-Nr, Name oder Teil davon'. Below it is a table with columns 'Name' and 'Pat.-Nr.'. The first row, 'Mustermann, Martina' with ID '123456', is highlighted in blue. A 'Neu einlesen' button with a refresh icon is located to the right of the search field. Below the table is a button labeled 'Diesen Patienten auswählen'. In the 'Untersuchungen' section, the header reads '123456 | Mustermann, Martina'. Below it is a table with columns 'Datum', 'Untersuchung', and 'Nr.'. The table contains five rows of examination records. At the bottom of the window, there are three buttons: 'Neuer Patient' and 'Patient Bearbeiten'. The status bar at the bottom shows 'Patientin/Patient Pat-Nr. - - -' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

Name	Pat.-Nr.
Mustermann, Martina	123456
Musterfrau, Mara	133578
Musterpatient, Paul	999781
Patient01, Patient1	1
Patient02, Patient2	2
Patient03, Patient3	3
Patient04, Patient4	4
Patient05, Patient5	5
Patient06, Patient6	6
Patient07, Patient7	7

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The screenshot displays the CMDfact 4 software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo and the 'dentaConcept' logo are visible. A toolbar contains various module icons: 'CMD CHECK', 'CMD brux', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. The 'CMD pain' icon is highlighted with a red circle and the number '1'. A hand icon is shown dragging the 'CMD pain' icon from the toolbar into a vertical timeline on the left. The timeline has a grey bar for the current day, circled with a red '2'. The 'CMD pain' icon is being placed into this bar, circled with a red '3'. A mouse cursor is shown clicking on the 'CMD pain' icon. The timeline is labeled 'Tag 1' with a calendar icon. The bottom of the interface shows 'Patientin/Patient Pat-Nr. ---' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot displays the CMDfact 4 software interface. At the top, there is a menu bar with 'Datei', 'Extras', and a question mark. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A horizontal toolbar contains several icons: a person icon, 'CMD CHECK', 'CMD brux', 'CMD status', 'CMD manu', 'CMD Pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. The 'CMD Pain' icon is highlighted with a blue circle and a red circle containing the number '2'. A hand icon is shown double-clicking this icon. Below the toolbar, a vertical sidebar contains four colored icons (grey, purple, red, yellow). The purple icon is highlighted with a blue circle and a red circle containing the number '1'. A hand icon is shown double-clicking this icon. To the right of the sidebar, a mouse icon is shown with a red circle containing the number '1' pointing to its button. Below the sidebar, a horizontal timeline arrow is labeled 'Tag 1' with a calendar icon. At the bottom of the window, the text 'Patientin/Patient Pat-Nr. ---' is on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' is on the right.

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The screenshot shows the CMDfact 4 software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A row of colorful icons represents different functions: CMD CHECK, CMD brux, CMD status, CMD manu, CMD pain, CMD tomo, CMD trace, CMD 3D, Diagnose Pilot, and CMD Arztbrief. The main area features a vertical timeline with five colored blocks (grey, purple, red, orange, yellow). A horizontal arrow points to the right from the bottom of the timeline, labeled 'Tag 1'. A calendar pop-up window is open, showing the month of July 2020. The date '12' is highlighted in blue. A red circle with the number '2' is positioned above the calendar, and a red circle with the number '1' is positioned below the 'Tag 1' label. The bottom status bar displays 'Patientin/Patient Pat-Nr. ---' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right.

Mo	Tu	We	Th	Fr	Sa	Su
		1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30	31		

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The screenshot displays the CMDfact 4 software interface. At the top, there is a menu bar with 'Datei', 'Extras', and a question mark. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A row of colorful icons represents different modules: a person icon, 'CMD CHECK', 'CMD brux', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. The main area shows a vertical timeline with five colored blocks (grey, purple, red, orange, yellow). A right-click context menu is open over the purple block, containing two options: 'Formblatt ausdrucken' and 'Untersuchung löschen'. A mouse cursor is positioned over the 'Untersuchung löschen' option. A horizontal arrow at the bottom of the timeline points to the right, with 'Tag 1' and a calendar icon at its start. The bottom status bar shows 'Patientin/Patient Pat-Nr. ---' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right.

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The screenshot shows the CMDfact 4 software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A row of icons represents various modules: CMD CHECK, CMD brux, CMD status, CMD manu, CMD pain, CMD tomo, CMD trace, CMD 3D, Diagnose Pilot, and CMD Arzt-brief. The main window displays a 'CMDstatus' window with a 'Palpationsbefunde' table. The table has columns for 'rechts' and 'links' and rows for various muscles and symptoms. A red circle '1' highlights the close button (X) in the window title bar. A red circle '2' highlights the 'CMDstatus' window title. A red arrow points from the 'CMDstatus' window to the 'Palpationsbefunde' table. The patient information at the bottom reads 'Mara Musterfrau Pat.-Nr. 133578' and 'Untersuchung vom 16.07.2017'. The footer includes 'Patientin/Patient Pat.-Nr. ---' and 'Praxis ... rium Hamburg-Eppendorf'.

	rechts	Schmerz	links
Missempfindung	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. masseter superfic.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. masseter prof.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Kiefergelenk lateral	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Kiefergelenk dorsal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. temporalis ant.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. temporalis med./post.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Suboccip.-/Nacken-M.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. trapezius	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. sternocleidomastoideus	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
infrahyoidale M.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
suprahyoidale M.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. pterygoideus med.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. digastricus venter post.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
M. pterygoideus lat.	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

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
The screenshot shows the CMDfact software interface. The top menu bar includes 'Datei', 'Extras', and '?'. The main header features the 'CMDfact' logo on the left and the 'dentaConcept' logo on the right. Below the header is a row of module icons: a person icon, 'CMD CHECK', 'CMD BRUX', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. Three red circles with numbers 1, 2, and 3 are placed below the 'CMD status', 'CMD pain', and 'Diagnose Pilot' icons, respectively, with lines pointing to them. The main content area is a large white box with the word 'Lieferumfang' (Scope of Delivery) in green text. A large grey arrow points to the right at the bottom of this box. The bottom status bar shows 'Patientin/Patient Pat-Nr. - - -' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right.

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The screenshot shows the CMDfact software interface. The menu bar includes 'Datei', 'Extras', and '?'. The main menu contains icons for 'CMD CHECK', 'CMD brux', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. Red circles with numbers 6 through 10 are placed over the 'CMD brux', 'CMD tomo', 'CMD trace', 'CMD 3D', and 'CMD Arzt-brief' icons respectively. The main content area displays 'in Entwicklung...' in blue text, with a large grey arrow pointing to the right. The status bar at the bottom shows 'Patientin/Patient Pat-Nr. - - -' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot shows the CMDstatus software interface. At the top, there are navigation icons for Stress, Ortho, Para, Palpation, Geräusche (highlighted), Mobilität, Okklusion1, Okklusion2, and Reaktion. Below these is a list of menu items under 'GELENKGERÄUSCHE', with 'Animation: Diskusverlagerung mit Reposition' selected. A text panel on the left provides a detailed description of disk displacement and its symptoms. On the right, a video player shows an anatomical illustration of the mandibular condyle and disk, with the title 'Kiefergelenk-Geräusche' and subtitle 'Diskusverlagerung mit Reposition'. A play button and volume control are visible below the video. On the far right, a vertical sidebar contains a switch for 'Anleitung' (instructions) and a 'Notizen' (notes) section.

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The screenshot displays the CMDstatus software interface. At the top, there are two window titles: 'CMDfact' and 'CMDstatus'. The main menu bar includes 'Datei', 'E', and 'CMD'. The 'CMDstatus' window has a header with the 'dentaConcept' logo and a toolbar with icons for 'Stress', 'Ortho', 'Para', 'Palpation', 'Geräusche', 'Mobilität', 'Okklusion1', 'Okklusion2', and 'Reaktion'. The 'Palpation' module is selected. Below the toolbar is a rich text editor with a 'Werkzeugleiste verbergen' button. The text area contains the heading '-- Notizen zu Palpationsbefunden' and the placeholder text 'Hier könnte eine Notiz stehen.' Below this, the text 'Sie können Sie formatieren | wie sie wollen.' is displayed in a bold, black serif font. A red circle with the number '2' points to the font size dropdown menu in the toolbar. On the right side of the window, there is a vertical sidebar with a 'Notizen' button (circled with a red '1') and an 'Anleitung' button. At the bottom of the window, the patient information 'Mara Musterfrau' and 'Pat-Nr. 133578' is shown, along with '232 Zeichen' and a checkbox for 'Notizen auf neue Seite drucken'. The footer of the application reads 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot displays the 'CMDfact' software interface. At the top, there's a menu bar with 'Datei' and 'Extras'. Below it, the 'CMDfact' logo is prominent. A navigation bar contains icons for 'CMD CHECK', 'CMD BRUX', 'CMD status', and 'CMD WA'. The main content area shows a detailed clinical report for 'Mara Mustermann', dated 27.02.2006. The report is titled 'Klinischer Funktionsstatus (4)' and includes sections for 'Stress', 'Orthopädisches Screening', 'Palpationsbefunde', 'Mobilität', 'Okklusionsbefunde 1', 'Okklusionsbefunde 2', 'Reaktionstests', and 'Notizen'. Each section contains various checkboxes, text boxes, and diagrams. For example, the 'Palpationsbefunde' section includes diagrams of the head and neck with red highlights indicating areas of concern. The 'Okklusionsbefunde 1' section features a grid for 'Öffnung UK' and 'Overbite (vert.)'. The 'Reaktionstests' section includes a 'Provokationstest (Klopp-Poulsen)' and a 'Resilienztest nach GERBER'. The 'Notizen' section at the bottom contains a box for additional notes. The interface is clean and professional, with a mix of text and visual aids.

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The screenshot shows the CMDfact software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A navigation bar contains several buttons: 'CMD CHECK' (highlighted with a red dashed box and a red arrow), 'CMD brux', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. The main content area displays a large orange splash screen for 'CMDCHECK Version 4 4.15'. The splash screen features a 3D anatomical illustration of a human head and neck with red dots indicating areas of interest. The 'dentaConcept' logo is visible at the bottom right of the splash screen. At the bottom of the software window, there is a status bar with 'Patientin/Patient Pat-Nr. ---' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right.

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The screenshot displays the CMDstatus software interface. At the top, there is a menu bar with 'Datei' and a question mark. Below the menu bar, the 'CMDstatus' logo is on the left, and the 'dentaConcept' logo is on the right. A horizontal toolbar contains nine icons representing different assessment categories: Stress (lightning bolt), Ortho (head with arrow), Para (knee joint), Palpation (head profile), Geräusche (ear), Mobilität (spine), Okklusion1 (teeth), Okklusion2 (teeth), and Reaktion (crossed arms). The 'Stress' icon is highlighted with a red circle labeled '1'. Below the toolbar, a large white panel is titled 'Stress' with a dark header and a white square icon. A red circle labeled '2' points to the text 'Übertrag aus dem Fragebogen Stressbelastung:'. Below this text is a text input field containing the placeholder 'hier eintragen'. To the right of the input field is the label 'Punkte'. Below the input field, the text 'Die Stressbelastung ist' is followed by three radio button options: 'extrem (> 200 Punkte)', 'erhöht (100-200 Punkte)', and 'normal (< 100 Punkte)'. At the bottom of the interface, patient information is displayed: 'Mara Musterfrau', 'Pat-Nr. 133578', and 'Untersuchung vom 16.07.2017'. The bottom status bar shows 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot shows the 'CMDstatus' software window. The title bar includes 'CMDfact' and 'CMDstatus'. The main menu bar has 'Datei' and '?'. The logo 'CMDstatus' and 'dentaConcept' are visible. A toolbar contains icons for Stress, Ortho, Para, Palpation, Geräusche, Mobilität, Okklusion1, Okklusion2, and Reaktion. The 'Ortho' icon is highlighted with a red circle labeled '1'. Below the toolbar, the 'Orthopädisches Screening' module is displayed, with its title bar also highlighted by a red circle labeled '2'. The module contains three main sections: 'Haltung seitlich' (Side Posture) with a diagram and measurement instructions (Maß: 2 Handdicken Abstand vom Lot), 'Haltung frontal' (Frontal Posture) with a diagram, and 'HWS' (Cervical Spine) with a diagram and three sub-sections: 'Inklination', 'Mittelstellung', and 'Reklination'. Each sub-section has three radio buttons labeled 'li', 'keine', and 're'. There are also checkboxes for 'Okklusionsänderung bei Haltungskorrektur' and 'Ligamentäre Schwäche'. The patient information at the bottom reads: 'Mara Musterfrau Pat-Nr. 133578' and 'Untersuchung vom 16.07.2017'. The footer includes 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot displays the CMDstatus software interface. At the top, there is a menu bar with 'Datei' and a question mark. Below the menu bar, the 'CMDstatus' logo is on the left, and the 'dentaConcept' logo is on the right. A row of icons represents different clinical areas: Stress (lightning bolt), Ortho (head with arrows), Para (knee joint), Palpation (head profile), Geräusche (ear), Mobilität (spine), Okklusion1 (teeth), Okklusion2 (teeth), and Reaktion (crossed arms). The 'Para' icon is highlighted with a red circle labeled '1'. Below the icons, a large panel titled 'Parafunktion' is shown, with a red circle labeled '2' pointing to its title bar. This panel contains a list of checkboxes for various conditions: Zungenindentationen, Wangenveränderungen, Abrasionen/Attritionen, keilförmige Defekte, Rezessionen, McCallsche Girlanden, Stillman-Spalten, and Zahnlockerungen... To the right of each checkbox is a text input field with the placeholder 'an den Zähnen...'. A red circle labeled '3' points to the first of these input fields. At the bottom of the interface, patient information is displayed: 'Mara Musterfrau', 'Pat-Nr. 133578', and 'Untersuchung vom 16.07.2017'. The bottom status bar includes 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot shows the 'CMDstatus' software interface. At the top, there are navigation icons for Stress, Ortho, Para, Palpation, Geräusche, Mobilität, Okklusion1, Okklusion2, and Reaktion. The 'Reaktionstests' window is open, displaying the following information:

- Reaktionstests** (with a close button 'X')
- Provokationstest KROUGH-POULSEN**
  - Radio buttons for 'pos.' (selected) and 'neg.'
  - Input field: 'Pressen 24,25' (with 'Ausführung' label)
- Resilienztest nach GERBER**
  - Two columns: 'rechts' and 'links', each with 'links' and 'rechts' sub-columns.
  - Each sub-column has a 'Shimstock fixiert bei mm Zinnfolie' label and an input field.
  - Radio buttons for 'Kompression', 'neutral', and 'Distraction' are present for each sub-column.
  - A checkbox with a red 'X' is checked, labeled 'Nicht anwendbar wegen Non-Okklusion'.

At the bottom of the window, the patient information is displayed: 'Mara Musterrau', 'Pat-Nr. 133578', and 'Untersuchung vom 16.07.2017'. The patient name is also visible in the footer of the main application window.

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The screenshot shows the CMDfact software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A horizontal toolbar contains several icons: a person icon, 'CMD CHECK', 'CMD brux', 'CMD status', 'CMD manu' (highlighted with a red circle and labeled '1'), 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. A red arrow points from the 'CMD manu' icon in the toolbar to a smaller 'CMD manu' icon in the main content area, which is also circled in red and labeled '2'. The main content area displays the 'CMDmanu' logo in a large, stylized font. At the bottom left of the main content area, there are three stacked icons: 'CMD status', 'CMD CHECK', and a date '16.07.2017'. At the bottom of the interface, there is a footer with the text 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right.

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The screenshot shows the CMDmanu 4 software interface. The top menu bar includes 'Datei', 'Ex', and 'CMDmanu'. Below the menu is a toolbar with icons for 'Öffnung', 'Schluss', 'Protrusion', 'Latero R', 'Latero L', 'Dyn.Kompr.', 'Stat.Kompr.', 'Translation', and 'Traktion'. A red dashed box highlights the 'Dyn.Kompr.' and 'Stat.Kompr.' icons, with a red circle containing the number '2' pointing to them. The main window displays 'Kompression in Ruhe (passive Kompression)'. Below this, there are two columns for 'Condylus rechts' and 'Condylus links'. Each column has checkboxes for 'Schmerz' and 'cranial'. Under 'Condylus rechts', 'dorso-cranio-lateral, dorso-cranial' and 'dorso-lateral, dorsal' are listed with checkboxes. Under 'Condylus links', 'cranial', 'dorso-cranio-lateral, dorso-cranial', and 'dorso-lateral, dorsal' are listed with checkboxes. The 'dorso-cranio-lateral, dorso-cranial' and 'dorso-lateral, dorsal' checkboxes are marked with a red 'X'. The patient information at the bottom reads: 'Mara Musterfrau Pat-Nr. 133578' and 'Untersuchung vom 16.07.2017'. The footer also includes 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot shows the CMDfact software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A horizontal toolbar contains several icons: a person icon, 'CMD CHECK', 'CMD DRUX', 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arztbrief'. The 'CMD pain' icon is circled with a red '1' and has a purple arrow pointing to it from the left. Below the toolbar, a large lightning bolt icon is centered. On the left side of the main content area, there is a vertical sidebar with four buttons: 'CMD pain', 'CMD manu', 'CMD status', and 'CMD CHECK'. The 'CMD pain' button is circled with a red '2' and has a purple arrow pointing to it from the top. At the bottom of the main content area, there is a date field showing '16.07.2017' and a large grey arrow pointing to the right. At the very bottom of the window, there is a footer with the text 'Patientin/Patient **Mara Musterfrau** Pat-Nr. 133578' on the left and 'Praxis **CMD-Centrum Hamburg-Eppendorf**' on the right.

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The screenshot shows the CMDpain 4 software interface. At the top, there are three main navigation buttons: 'Hauptbeschwerden' (Main Complaints), 'Schmerzen' (Pain), and 'Verlauf' (Course). The 'Hauptbeschwerden' section is active, displaying a text input field with the text 'Schmerzen KG li'. Below this is a pain scale section titled 'Schmerzstärke insgesamt' (Overall Pain Intensity). The scale ranges from '0 = gar nicht' (0 = not at all) to 'unerträglich = 10' (unbearable = 10). A purple bar indicates the current pain level, which is set to 8. A hand icon is shown interacting with the scale, pointing to the value 8. At the bottom of the interface, the patient's name 'Mara Musterfrau' and patient number 'Pat-Nr. 133578' are displayed, along with the date '16.07.2017' and the practice name 'Praxis CMD-Centrum Hamburg-Eppendorf'. A checkbox labeled 'Werte anzeigen' (Show values) is checked. The interface also includes a sidebar with 'Notizen' (Notes) and 'Anleitung' (Instructions) sections.

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The screenshot shows the CMDpain software interface. The main window displays the patient's name, 'Mara Musterfrau', and patient number, 'Pat-Nr. 133578'. The assessment date is '16.07.2017'. The interface is divided into three main sections: 'Hauptbeschwerden', 'Schmerzen', and 'Verlauf'. The 'Schmerzen' section is active and shows four categories of pain, each with a horizontal bar chart and a numerical value. A red circle highlights the value '7' for 'Schmerzen an Schläfen oder Wangen'. The other categories and their values are: 'Schmerzen an den Zähnen' (4), 'Schmerzen im Kiefergelenk oder am Ohr' (6), and 'Schmerzen in Nacken oder Rücken' (3). The interface also includes a 'Notizen' sidebar on the right and a 'Werte anzeigen' checkbox at the bottom.

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The screenshot displays the CMDpain 4 software interface. At the top, there are three main navigation buttons: 'Hauptbeschwerden', 'Schmerzen', and 'Verlauf'. The 'Verlauf' button is currently selected. Below these buttons is a large graph titled 'Verlauf von Hauptbeschwerden und lokalen Schmerzen'. The y-axis is labeled 'Schmerzen' and ranges from 0 to 10. The x-axis is labeled 'Untersuchungstermine'. A red circle with the number '1' is placed over the y-axis at the value 8. Below the graph, there are five checkboxes for different pain categories: 'Gesamtschmerz', 'Schläfen/Wangen', 'Zähne', 'Kiefergelenke', and 'Nacken/Rücken'. A red circle with the number '2' is placed over the 'Schläfen/Wangen' checkbox. At the bottom of the interface, patient information is displayed: 'Mara Musterfrau', 'Pat-Nr. 133578', and 'Untersuchung vom 16.07.2017'. The footer of the window shows 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot displays the CMDfact software interface. At the top, there is a menu bar with 'Datei', 'Extras', and '?'. Below the menu bar, the 'CMDfact' logo is on the left and the 'dentaConcept' logo is on the right. A horizontal toolbar contains several icons: a person icon, 'CMD CHECK', 'CMD briix', 'CMD status', 'CMD manu', 'CMD pain' (circled with a red '1'), 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot' (highlighted in green), and 'CMD Arzt-brief'. A green arrow points from the 'Diagnose Pilot' icon in the toolbar to a vertical sidebar on the left. This sidebar contains a list of icons: 'Diagnose Pilot' (green), 'CMD pain' (purple), 'CMD manu' (red), 'CMD status' (red), and 'CMD CHECK' (yellow). A red circle with the number '2' is around the 'Diagnose Pilot' icon in the sidebar. A large question mark is centered in the main workspace area. At the bottom of the interface, a date field shows '16.07.2017'. The footer contains the text 'Patientin/Patient **Mara Musterfrau** Pat-Nr. 133578' on the left and 'Praxis **CMD-Centrum Hamburg-Eppendorf**' on the right.

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The screenshot displays the 'CMD DiagnosePilot' software interface. The top navigation bar includes buttons for '(Initial-) Diagnosen', 'Neben- + Differential-Diagnosen', 'Diagnose-Cinema', 'Maßnahmen, Konsile', and 'Export'. The main content area is divided into two columns. The left column, titled 'Initialdiagnosen', lists various conditions with checkboxes for selection and columns for 'R' and 'L' side findings. The right column, titled 'Anamnesen und Befunde', is further divided into 'Leitsymptome', 'passende Befunde', and 'NICHT passende Befunde'. A sidebar on the left contains a vertical menu with icons for 'DiagnosePilot', 'CMD', 'Maßnahmen', 'Status', and 'CMD'. The patient information at the bottom identifies 'Michaela Mustermann' (Pat-Nr. 0815) and 'Mara Musterfrau' (Pat-Nr. 133578). The status bar shows 'CMDstatus 12.5.2017 | DiagnosePilot 12.05.2017' and the practice name 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot shows the CMD DiagnosePilot software interface. At the top, there are menu options 'Datei' and 'Extras'. Below the title bar, the text 'CMD DiagnosePilot' and the 'dentaConcept' logo are visible. A navigation bar contains five green buttons: '(Initial-) Diagnosen', 'Neben- + Differential-Diagnosen', 'Diagnose-Cinema', 'Maßnahmen, Konsile', and 'Export, Drucken'. A red circle with the number '1' highlights the '(Initial-) Diagnosen' button. The main content area features a large white box with the heading 'Kurzanleitung zum CMD DiagnosePilot'. Below the heading, a paragraph states: 'Der CMD DiagnosePilot fasst die erhobenen Befunde zusammen und erleichtert Ihnen die Auswahl der zutreffenden Diagnosen. CMDfact und der CMD DiagnosePilot sind keine Medizinprodukte, alle Entscheidungen über Diagnosen treffen allein Sie!'. A red circle with the number '2' points to the bottom of this text box. To the right of the text is a laptop displaying the software interface. A red circle with the number '3' points to the right side of the laptop. At the bottom of the laptop screen, there is a checkbox labeled 'Beim Start anzeigen' and a 'Schließen' button. The bottom status bar of the software shows 'Patientin/Patient Mara Musterfrau Pat-Nr. 133578' and 'Praxis CMD-Centrum Hamburg-Eppendorf'. A date '16.07.20' is visible in the bottom left corner of the software window.

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The screenshot displays the CMDfact software interface. The main window has a menu bar with 'Datei', 'Extras', and '?'. The 'Extras' menu is open, showing options like 'Praxisdaten / Aktivierungscodes', 'Datenbank-Verbindung', 'VDDS anmelden', 'Einstellungen', 'Services', and 'Import von CMDfact3'. Below the menu is a toolbar with various module icons: 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arzt-brief'. A red circle with the number '2' highlights the 'CMD manu' icon. In the center, a dialog box titled 'Freischaltung' is open, displaying 'FREIGESCHALTET' with a red circle and the number '1' next to it. Below this, it says 'Aktivierungscode und Registrier-Name waren korrekt' and has an 'OK' button. At the bottom of the main window, there is a status bar with 'Patientin/Patient - - - Pat-Nr. - - -' on the left and 'Praxis CMD-Centrum Hamburg-Eppendorf' on the right. A large grey arrow points from the dialog box towards the right side of the screen.

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The screenshot shows the CMDfact software interface. The 'Extras' menu is open, and the 'VDDS anmelden' option is highlighted with a red circle and the number '1'. A red arrow points from this menu item to the 'CMD VDDS' button in the top navigation bar. Below this, a dialog box titled 'CMDfact : VDDS aktivieren' is displayed. The dialog box has a black header with the text 'VDDSmedia-Schnittstelle'. The main content area contains the following text:

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At the bottom of the dialog box, there are two buttons: 'Schließen' and 'In VDDS anmelden'. The 'In VDDS anmelden' button is highlighted with a red circle and the number '2'. The status bar at the bottom of the main window shows 'Patientin/Patient - - - Pat-Nr. - - -' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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The screenshot displays the CMDfact 4 software interface. The main window has a menu bar with 'Datei', 'Extras', and '?'. The 'Extras' menu is open, showing options like 'Praxisdaten / Aktivierungscodes', 'Datenbank-Verbindung', 'VDDS anmelden', 'Einstellungen', 'Services', and 'Import von CMDfact3'. The 'Einstellungen' dialog box is open, showing settings for 'Dateiname der Ausdrucke' and 'Farbe der Messwerte'. The 'Dateiname der Ausdrucke' section includes input fields for 'Patienten-Nummer', 'Vorname des/der Patienten/in', 'Nachname des/der Patienten/in', 'Datum der Untersuchung', and 'Name der Untersuchung/des Moduls'. The 'Muster' field contains the text '12345 Berger Eva CMDstatus 23-23-2003.pdf', with a red circle '1' highlighting it and red arrows pointing to the individual components of the filename. Below the 'Muster' field are buttons for 'Neu beginnen', 'Standard', and 'zurücksetzen'. A red circle '2' highlights the 'Standard' button. The 'Farbe der Messwerte' section includes a color selection dropdown set to '#4d66cc' and a 'zurücksetzen' button. At the bottom of the dialog are 'Abbruch' and 'Sichern' buttons. The background shows the main software interface with the 'dentaConcept' logo and various module icons like 'CMD status', 'CMD manu', 'CMD pain', 'CMD tomo', 'CMD trace', 'CMD 3D', 'Diagnose Pilot', and 'CMD Arztbrief'. The status bar at the bottom shows 'Patientin/Patient - - - Pat-Nr. - - -' and 'Praxis CMD-Centrum Hamburg-Eppendorf'.

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